

Elastase 1 in serum

Template:Zkontrolováno thumb|Markery akutní pankreatitidy

Introduction

- **Elastase-1** (EC 3.4.21.11) is secreted as proelastase and activated by trypsinem.
- We distinguish:
 1. elastase-1 (Mr 30 000; anodic fraction), which occurs in serum in free form and in complex with α_1 -proteinase inhibitor
 2. elastase-2 (Mr 25 000; cathodic fraction).
- The enzyme is secreted by the pancreatic juice into the duodena and the protein sequence selected for [immunochemical detection]] při stanovení elastase-1 in stool.
- In inflammatory processes, there is also a retrograde release into the bloodstream and quantification of serum levels of human pancreatic elastase is a suitable marker of acute pancreatitis, such as post-ERCP, and pancreatic cancer.

Characteristics

- Serum elastase-1 (determined by ELISA) has the highest of all enzymes for pancreatic cancer specificity and sensitivity.
- Serum elastase-1 levels are elevated in both acute and chronic recurrent pancreatitis, and the increase persists longer and correlates better with clinical status than α -amylase.
- The reference range (for RIA methodology) is 1.3–4.3 $\mu\text{g/l}$. The latest study in the determination of serum elastase 1 is based on a latex immunoassay.

Links

Source

- with the permission of the author taken from KOCNA, Petr. *GastroLab: MiniEncyclopedia of laboratory methods in gastroenterology* [online]. © 2002. Last revision 2011-01-08, [cit. 2011-03-04]. <<http://www1.lf1.cuni.cz/~kocna/qlab/glency1.htm> >.

References

- KAPETANOS, D, et al. The value of serum amylase and elastase measurements in the prediction of post-ERCP acute pancreatitis. *Hepatogastroenterology*. 2007, vol 54, no. 74, pp. 556-60, ISSN 0172-6390. PMID: 17523321 .
- WU, D, et al. Evaluating the role of serum elastase 1 in the diagnosis of pancreatic cancer. *Chin J Dig Dis*. 2006, vol 7, no. 2, pp. 117-20, ISSN 1443-9611 (Print), 1443-9573 (Electronic). PMID: 16643340 .
- WILSON, RB, et al. Serum elastase in the diagnosis of acute pancreatitis: a prospective study. *ANZ J Surg*. 2005, vol 75, no. 3, pp. 152-6, ISSN 1445-1433 (Print), 1445-2197 (Electronic). PMID: 15777396 .
- KEIM, V, et al. Evaluation of Pankrin, a new serum test for diagnosis of acute pancreatitis. *Clin Chim Acta*. 2003, vol 332, no. 1-2, pp. 45-50, ISSN 0009-8981 (Print), 1873-3492 (Electronic). PMID: 12763279 .
- KATSANOS, KH, et al. The true value of serum elastase-1 in endoscopic retrograde cholangiopancreatography (ERCP). *Eur J Intern Med*. 2002, vol 13, no. 5, pp. 329-335, ISSN 0953-6205 (Print), 1879-0828 (Electronic). PMID: 12144913 .
- OHDE, K, et al. Determination of elastase 1 in plasma or serum by latex turbidometric immunoassay with automatic analyzer. *Rinsho Byori*. 2002, vol 50, no. 4, pp. 376-80, ISSN 0047-1860. PMID: 12014017 .
- MILLSON, CE, et al. A prospective study of serum pancreatic elastase-1 in the diagnosis and assessment of acute pancreatitis. *Scand J Gastroenterol*. 1998, vol. 33, no. 6, pp. 664-8, ISSN 0036-5521 (Print), 1502-7708 (Electronic). PMID: 9669641 .
- IKEI, S, et al. Blood concentrations of polymorphonuclear leucocyte elastase and interleukin-6 are indicators for the occurrence of multiple organ failures at the early stage of acute pancreatitis. *J Gastroenterol Hepatol*. 1998, vol 13, no. 12, pp. 1274-83, ISSN 0815-9319 (Print), 1440-1746 (Electronic). PMID: 9918438 .
- KEIM, V, et al. Polyclonal pancreatic elastase assay is superior to monoclonal assay for diagnosis of acute pancreatitis. *Clin Chem*. 1997, vol 43, no. 12, pp. 2339-44, ISSN 0009-9147 (Print), 1530-8561 (Electronic). PMID: 9439452 .

- IWASE, K, et al. Serial changes in plasma levels of pancreatic elastase 1 after pancreatic surgeries. *Res Exp Med (Berl)*. 1995, vol 195, no. 2, pp. 93-100, ISSN 0300-9130. PMID: 7659839 .
- SATAKE, K, et al. CA19-9 as a screening and diagnostic tool in symptomatic patients: the Japanese experience. *Pancreas*. 1994, vol. 9, no. 6, pp. 703-6, ISSN 0885-3177 (Print), 1536-4828 (Electronic). PMID: 7846012 .
- HOMMA, T, et al. The study of the mass screening of persons without symptoms and of the screening of outpatients with gastrointestinal complaints or icterus for pancreatic cancer in Japan, using CA19-9 and elastase-1 or ultrasonography. *Int J Pancreatol*. 1991, vol. 9, pp. 119-24, ISSN 0169-4197. PMID: 1744437 .
- LINDER, R, et al. Pancreatic elastase 1 after pancreatic transplantation. *Pancreas*. 1991, vol. 6, no. 1, pp. 31-6, ISSN 0885-3177 (Print), 1536-4828 (Electronic). PMID: 1704632 .
- PEZZILLI, R, et al. Serum pancreatic enzymes in the diagnosis of carcinoma of the pancreas. *G Ital Oncol*. 1989, vol. 9, no. 1, pp. 7-13, ISSN 0392-128X. PMID: 2468605 .